

产品名称: ALS2CR11 Rabbit Polyclonal Antibody
产品货号: APRab06811



产品概述 (Summary)

产品名称 (Production Name)	ALS2CR11 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ALS2CR11
别名 (Alternative Names)	ALS2CR11; Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 11 protein
基因 ID (Gene ID)	151254.0
蛋白 ID (SwissProt ID)	Q53TS8.The antiserum was produced against synthesized peptide derived from human ALS2CR11. AA range:491-540

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	72kDa

产品名称: ALS2CR11 Rabbit Polyclonal Antibody
产品货号: APRab06811




研究背景 (Background)


An autosomal recessive form of juvenile amyotrophic lateral sclerosis was originally mapped to a region of chromosome 2 that includes this gene. The encoded protein contains a calcium-dependent membrane targeting C2 domain. This domain is often found in proteins that are involved in membrane trafficking and signal transduction. [provided by RefSeq, Jun 2016],

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of ALS2CR11 Antibody. The lane on the right is blocked with the ALS2CR11 peptide.



Western blot analysis of the lysates from COLO205 cells using ALS2CR11 antibody.

产品名称: ALS2CR11 Rabbit Polyclonal Antibody
产品货号: APRab06811



Western Blot analysis of 3T3 cells using ALS2CR11 Polyclonal Antibody

注意事项 (Note)

For research use only .