

产品名称: ALS2CR13 Rabbit Polyclonal Antibody
产品货号: APRab06812



产品概述 (Summary)

产品名称 (Production Name)	ALS2CR13 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FAM117B
别名 (Alternative Names)	FAM117B; ALS2CR13; Protein FAM117B; Amyotrophic lateral sclerosis 2 chromosomal region candidate gene 13 protein
基因 ID (Gene ID)	150864.0
蛋白 ID (SwissProt ID)	Q6P1L5.The antiserum was produced against synthesized peptide derived from human ALS2CR13. AA range:411-460

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:10000
蛋白分子量 (Molecular Weight)	60kDa

产品名称: ALS2CR13 Rabbit Polyclonal Antibody
产品货号: APRab06812




研究背景 (Background)


miscellaneous:ALS2CR13 is mapped in the genomic region covering the complete candidate region for Amyotrophic lateral sclerosis 2 (ALS2).,sequence caution:Contaminating sequence. Sequence of unknown origin in the N-terminal part.,miscellaneous:ALS2CR13 is mapped in the genomic region covering the complete candidate region for Amyotrophic lateral sclerosis 2 (ALS2).,sequence caution:Contaminating sequence. Sequence of unknown origin in the N-terminal part.,

研究领域 (Research Area)

图片 (Image Data)



Western Blot analysis of HELA cells using ALS2CR13 Polyclonal Antibody diluted at 1: 500



Western Blot analysis of HELA cells using ALS2CR13 Polyclonal Antibody diluted at 1: 500

注意事项 (Note)

For research use only .