

产品名称: ALS2CR8 Rabbit Polyclonal Antibody
产品货号: APRab06815



产品概述 (Summary)

产品名称 (Production Name)	ALS2CR8 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ALS2CR8 ALS2CR8; CARF; Amyotrophic lateral sclerosis 2 chromosomal region
别名 (Alternative Names)	candidate gene 8 protein; Calcium-response factor; CaRF; Testis development protein NYD-SP24
基因 ID (Gene ID)	79800.0
蛋白 ID (SwissProt ID)	Q8N187.The antiserum was produced against synthesized peptide derived from human ALS2CR8. AA range:311-360

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	85kDa

产品名称: ALS2CR8 Rabbit Polyclonal Antibody
产品货号: APRab06815




研究背景 (Background)

function:May be a transcription factor.,function:May be a transcription factor.,


研究领域 (Research Area)

Cell Biology

图片 (Image Data)




Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue, using ALS2CR8 Antibody. The picture on the right is blocked with the synthesized peptide.



Western blot analysis of lysates from HUVEC cells, using ALS2CR8 Antibody. The lane on the right is blocked with the synthesized peptide.

产品名称: ALS2CR8 Rabbit Polyclonal Antibody
产品货号: APRab06815



Western Blot analysis of various cells using ALS2CR8 Polyclonal Antibody cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventbiotech, MN, USA) .

注意事项 (Note)

For research use only .