

产品名称: AMBRA1 Rabbit Polyclonal Antibody
产品货号: APRab06821



产品概述 (Summary)

产品名称 (Production Name)	AMBRA1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	AMBRA1 KIAA1736
别名 (Alternative Names)	autophagy/beclin-1 regulator 1
基因 ID (Gene ID)	55626.0
蛋白 ID (SwissProt ID)	Q9C0C7.Synthetic Peptide of AMBRA1

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	140kDa

研究背景 (Background)

function:Regulates autophagy and development of the nervous system. Involved in autophagy in controlling protein

产品名称: AMBRA1 Rabbit Polyclonal Antibody
产品货号: APRab06821




turnover during neuronal development, and in regulating normal cell survival and proliferation.,similarity:Contains 3 WD repeats.,subunit:Interacts with BECN1. Probably forms a complex with BECN1 and PIK3C3.,function:Regulates autophagy and development of the nervous system. Involved in autophagy in controlling protein turnover during neuronal development, and in regulating normal cell survival and proliferation.,similarity:Contains 3 WD repeats.,subunit:Interacts with BECN1. Probably forms a complex with BECN1 and PIK3C3.,

研究领域 (Research Area)

Autophagy

图片 (Image Data)



Western Blot analysis of 293T 3T3 cells using AMBRA1 Polyclonal Antibody diluted at 1:500. Secondary antibody was diluted at 1:20000

注意事项 (Note)

For research use only .