

**产品名称: Annexin VI Rabbit Polyclonal Antibody**  
**产品货号: AP Rab06927**



## 产品概述 (Summary)

|                        |                                       |
|------------------------|---------------------------------------|
| 产品名称 (Production Name) | Annexin VI Rabbit Polyclonal Antibody |
| 描述 (Description)       | Rabbit polyclonal Antibody            |
| 宿主 (Host)              | Rabbit                                |
| 应用 (Application)       | WB,IHC,ICC/IF,ELISA                   |
| 种属反应性 (Reactivity)     | Human,Rat,Mouse                       |

## 产品性能 (Performance)

|                     |                                                                                                   |
|---------------------|---------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                      |
| 修饰 (Modification)   | Unmodified                                                                                        |
| 同种型 (Isotype)       | IgG                                                                                               |
| 克隆 (Clonality)      | Polyclonal                                                                                        |
| 形式 (Form)           | Liquid                                                                                            |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.          |
| 储存溶液 (Buffer)       | Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N. |
| 纯化方式 (Purification) | Affinity purification                                                                             |

## 免疫原信息 (Immunogen)

|                        |                                                                                                                                                    |
|------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | ANXA6                                                                                                                                              |
| 别名 (Alternative Names) | ANXA6; ANX6; Annexin A6; 67 kDa calelectrin; Annexin VI; Annexin-6; Calphobindin-II; CPB-II; Chromobindin-20; Lipocortin VI; Protein III; p68; p70 |
| 基因 ID (Gene ID)        | 309.0                                                                                                                                              |
| 蛋白 ID (SwissProt ID)   | P08133.The antiserum was produced against synthesized peptide derived from human Annexin A6. AA range:1-50                                         |

## 产品应用 (Application)

|                          |                                                                          |
|--------------------------|--------------------------------------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:100-1:300,ELISA 1:10000-1:20000 |
| 蛋白分子量 (Molecular Weight) | 75kDa                                                                    |

**产品名称: Annexin VI Rabbit Polyclonal Antibody**  
**产品货号: AP Rab06927**




## 研究背景 (Background)

Annexin VI belongs to a family of calcium-dependent membrane and phospholipid binding proteins. Several members of the annexin family have been implicated in membrane-related events along exocytotic and endocytotic pathways. The annexin VI gene is approximately 60 kbp long and contains 26 exons. It encodes a protein of about 68 kDa that consists of eight 68-amino acid repeats separated by linking sequences of variable lengths. It is highly similar to human annexins I and II sequences, each of which contain four such repeats. Annexin VI has been implicated in mediating the endosome aggregation and vesicle fusion in secreting epithelia during exocytosis. Alternatively spliced transcript variants have been described. [provided by RefSeq, Aug 2010], caution: The sequence shown here is derived from an Ensembl automatic analysis pipeline and should be considered as preliminary data., domain: A pair of annexin repeats may form one binding site for calcium and phospholipid., function: May associate with CD21. May regulate the release of Ca(2+) from intracellular stores., induction: By EBV., miscellaneous: Seems to bind one calcium ion with high affinity., PTM: Phosphorylated in response to growth factor stimulation., similarity: Belongs to the annexin family., similarity: Contains 4 annexin repeats., similarity: Contains 8 annexin repeats., subcellular location: Identified by mass spectrometry in melanosome fractions from stage I to stage IV.,

## 研究领域 (Research Area)


Cell Biology

## 图片 (Image Data)




Immunofluorescence analysis of HeLa cells, using Annexin A6 Antibody. The picture on the right is blocked with the synthesized peptide.


产品名称: Annexin VI Rabbit Polyclonal Antibody  
产品货号: APRab06927



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using Annexin A6 Antibody. The picture on the right is blocked with the synthesized peptide.




Western blot analysis of lysates from 293 cells, using Annexin A6 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using Annexin VI Polyclonal Antibody diluted at 1: 2000

产品名称: Annexin VI Rabbit Polyclonal Antibody  
产品货号: APRab06927



Western Blot analysis of 293 cells using Annexin VI Polyclonal Antibody diluted at 1: 2000

#### 注意事项 (Note)

For research use only .