

产品名称: Annexin X Rabbit Polyclonal Antibody
产品货号: APRab06929



产品概述 (Summary)

产品名称 (Production Name)	Annexin X Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ANXA10
别名 (Alternative Names)	ANXA10; ANX14; Annexin A10; Annexin-10; Annexin-14
基因 ID (Gene ID)	11199.0
蛋白 ID (SwissProt ID)	Q9UJ72.The antiserum was produced against synthesized peptide derived from human ANXA10. AA range:261-310

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	38kDa

研究背景 (Background)


产品名称: Annexin X Rabbit Polyclonal Antibody
产品货号: AP Rab06929




This gene encodes a member of the annexin family. Members of this calcium-dependent phospholipid-binding protein family play a role in the regulation of cellular growth and in signal transduction pathways. The function of this gene has not yet been determined. [provided by RefSeq, Jul 2008],similarity:Belongs to the annexin family.,similarity:Contains 4 annexin repeats.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of the lysates from Jurkat cells using ANXA10 antibody.



Western Blot analysis of various cells using Annexin X Polyclonal Antibody

注意事项 (Note)

For research use only .