

产品名称: Anoctamin-9 Rabbit Polyclonal Antibody
产品货号: APRab06935



产品概述 (Summary)

产品名称 (Production Name)	Anoctamin-9 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ANO9
别名 (Alternative Names)	ANO9; PIG5; TMEM16J; TP53I5; Anoctamin-9; Transmembrane protein 16J; Tumor protein p53-inducible protein 5; p53-induced gene 5 protein
基因 ID (Gene ID)	338440.0
蛋白 ID (SwissProt ID)	A1A5B4. The antiserum was produced against synthesized peptide derived from human TM16J. AA range: 481-530

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:300, ICC/IF 1:200-1:1000, ELISA 1:5000-1:10000
蛋白分子量 (Molecular Weight)	

产品名称: Anoctamin-9 Rabbit Polyclonal Antibody
产品货号: APRab06935




研究背景 (Background)


function:May act as a calcium-activated chloride channel.,sequence caution:Intron retention.,similarity:Belongs to the anoctamin family.,function:May act as a calcium-activated chloride channel.,sequence caution:Intron retention.,similarity:Belongs to the anoctamin family.,

研究领域 (Research Area)

图片 (Image Data)



Immunofluorescence analysis of HepG2 cells, using TM16J Antibody. The picture on the right is blocked with the synthesized peptide.



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using TM16J Antibody. The picture on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .