

产品名称: ApoF Rabbit Polyclonal Antibody
产品货号: APRab07040



产品概述 (Summary)

产品名称 (Production Name)	ApoF Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	APOF
别名 (Alternative Names)	APOF; Apolipoprotein F; Apo-F; Lipid transfer inhibitor protein; LTIP
基因 ID (Gene ID)	319.0
蛋白 ID (SwissProt ID)	Q13790.The antiserum was produced against synthesized peptide derived from human APOF. AA range:221-270

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	38kDa

研究背景 (Background)


产品名称: ApoF Rabbit Polyclonal Antibody
产品货号: APRab07040




The product of this gene is one of the minor apolipoproteins found in plasma. This protein forms complexes with lipoproteins and may be involved in transport and/or esterification of cholesterol. [provided by RefSeq, Jul 2008],function:Minor apolipoprotein that associates with LDL. Inhibits cholesteryl ester transfer protein (CETP) activity and appears to be an important regulator of cholesterol transport. Also associates to a lesser degree with VLDL, Apo-AI and Apo-AII.,tissue specificity:Liver. Secreted in plasma.,

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of lysates from HeLa, Jurkat, COLO, and A549 cells, using APOF Antibody. The lane on the right is blocked with the synthesized peptide.




Western blot analysis of the lysates from COLO205 cells using APOF antibody.

产品名称: ApoF Rabbit Polyclonal Antibody
产品货号: APRab07040



Western Blot analysis of various cells using ApoF Polyclonal Antibody diluted at 1: 1000



Western Blot analysis of A549 cells using ApoF Polyclonal Antibody diluted at 1: 1000

注意事项 (Note)

For research use only .