

产品名称: ApoOL Rabbit Polyclonal Antibody
产品货号: APRab07054



产品概述 (Summary)

产品名称 (Production Name)	ApoOL Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	APOOL
别名 (Alternative Names)	APOOL; CXorf33; FAM121A; Apolipoprotein O-like; Protein FAM121A
基因 ID (Gene ID)	139322.0
蛋白 ID (SwissProt ID)	Q6UXV4.The antiserum was produced against synthesized peptide derived from human APOOL. AA range:219-268

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	30kDa

研究背景 (Background)


产品名称: ApoOL Rabbit Polyclonal Antibody
产品货号: APRab07054




This gene encodes a protein which contains an apolipoprotein O superfamily domain. This domain is found on proteins in circulating lipoprotein complexes. [provided by RefSeq, Sep 2011],similarity:Belongs to the apolipoprotein O family.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from COLO and NIH/3T3 cells, using APOOL Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of SH-SY5Y 293T 3T3 lysis using ApoOL antibody. Antibody was diluted at 1:1000

注意事项 (Note)

For research use only .