

产品名称: Arfaptin 1 Rabbit Polyclonal Antibody
产品货号: APRab07104



产品概述 (Summary)

产品名称 (Production Name)	Arfaptin 1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ARFIP1
别名 (Alternative Names)	ARFIP1; Arfaptin-1; ADP-ribosylation factor-interacting protein 1
基因 ID (Gene ID)	27236.0
蛋白 ID (SwissProt ID)	P53367.The antiserum was produced against synthesized peptide derived from human ARFIP1. AA range:271-320

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	40kDa

研究背景 (Background)


产品名称: Arfaptin 1 Rabbit Polyclonal Antibody
产品货号: APRab07104




function: Putative target protein of ADP-ribosylation factor., similarity: Contains 1 AH domain., subunit: Interacts with non-myristoylated GTP-bound ARF3, but not to GDP-bound ARF3, and also to ARF1. Binds with lower affinity to ARF5 and with very little affinity to ARF6., tissue specificity: Ubiquitously expressed. Higher levels in liver, pancreas, placenta, skeletal muscle and heart., function: Putative target protein of ADP-ribosylation factor., similarity: Contains 1 AH domain., subunit: Interacts with non-myristoylated GTP-bound ARF3, but not to GDP-bound ARF3, and also to ARF1. Binds with lower affinity to ARF5 and with very little affinity to ARF6., tissue specificity: Ubiquitously expressed. Higher levels in liver, pancreas, placenta, skeletal muscle and heart.,

研究领域 (Research Area)

图片 (Image Data)




Immunofluorescence analysis of A549 cells, using ARFIP1 Antibody. The picture on the right is blocked with the synthesized



Western blot analysis of lysates from Jurkat cells, using ARFIP1 Antibody. The lane on the right is blocked with the synthesized peptide.

产品名称: Arfaptin 1 Rabbit Polyclonal Antibody
产品货号: APRab07104



Immunohistochemical analysis of paraffin-embedded Human thymus gland. Antibody was diluted at 1:100 (4°,overnight) . High-pressure and temperature Tris-EDTA,pH8.0 was used for antigen retrieval. Negative contrl (right) obtaned from antibody was pre-absorbed by immunogen peptide.

注意事项 (Note)

For research use only .