

产品名称: ARHGEF19 Rabbit Polyclonal Antibody
产品货号: APRab07129



产品概述 (Summary)

产品名称 (Production Name)	ARHGEF19 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ARHGEF19
别名 (Alternative Names)	ARHGEF19; Rho guanine nucleotide exchange factor 19; Ephexin-2
基因 ID (Gene ID)	128272.0
蛋白 ID (SwissProt ID)	Q8IW93.The antiserum was produced against synthesized peptide derived from human ARHGEF19. AA range:251-300

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	90kDa

研究背景 (Background)

产品名称: ARHGEF19 Rabbit Polyclonal Antibody
产品货号: AP Rab07129




Guanine nucleotide exchange factors (GEFs) such as ARHGEF19 accelerate the GTPase activity of Rho GTPases (see RHOA, MIM 165390).[supplied by OMIM, Dec 2008],function:Acts as guanine nucleotide exchange factor (GEF) for RhoA GTPase.,similarity:Contains 1 DH (DBL-homology) domain.,similarity:Contains 1 PH domain.,similarity:Contains 1 SH3 domain.,


研究领域 (Research Area)

Regulation of Actin Dynamics; AMPK

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma, using ARHGEF19 Antibody. The picture on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using ARHGEF19 Polyclonal Antibody

注意事项 (Note)

For research use only .