

产品名称: ARMC6 Rabbit Polyclonal Antibody
产品货号: APRab07147



产品概述 (Summary)

产品名称 (Production Name)	ARMC6 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, ELISA
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ARMC6
别名 (Alternative Names)	ARMC6; Armadillo repeat-containing protein 6
基因 ID (Gene ID)	93436.0
蛋白 ID (SwissProt ID)	Q6NXX6. The antiserum was produced against synthesized peptide derived from human ARMC6. AA range: 431-480

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:200-1:1000, ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	54kDa

研究背景 (Background)


产品名称: ARMC6 Rabbit Polyclonal Antibody
产品货号: APRab07147




The function of this gene's protein product has not been determined. A related protein in mouse suggests that this protein has a conserved function. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Nov 2010],PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Belongs to the ARMC6 family.,similarity:Contains 4 ARM repeats.,

研究领域 (Research Area)

图片 (Image Data)




Immunofluorescence analysis of A549 cells, using ARMC6 Antibody. The picture on the right is blocked with the synthesized



Western blot analysis of lysates from HeLa and COLO cells, using ARMC6 Antibody. The lane on the right is blocked with the synthesized peptide.

产品名称: ARMC6 Rabbit Polyclonal Antibody
产品货号: APRab07147



Western Blot analysis of various cells using ARMC6 Polyclonal Antibody

注意事项 (Note)

For research use only .