

产品名称: ARMCX1 Rabbit Polyclonal Antibody
产品货号: APRab07149



产品概述 (Summary)

产品名称 (Production Name)	ARMCX1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ARMCX1
别名 (Alternative Names)	ARMCX1; ALEX1; AD032; Armadillo repeat-containing X-linked protein 1; ARM protein lost in epithelial cancers on chromosome X 1; Protein ALEX1
基因 ID (Gene ID)	51309.0
蛋白 ID (SwissProt ID)	Q9P291.The antiserum was produced against synthesized peptide derived from human ARMX1. AA range:1-50

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	49kDa

产品名称: ARMCX1 Rabbit Polyclonal Antibody
产品货号: APRab07149




研究背景 (Background)


armadillo repeat containing, X-linked 1(ARMCX1) Homo sapiens This gene encodes a member of the ALEX family of proteins and may play a role in tumor suppression. The encoded protein contains a potential N-terminal transmembrane domain and two Armadillo (arm) repeats. Other proteins containing the arm repeat are involved in development, maintenance of tissue integrity, and tumorigenesis. This gene is closely localized with other family members, including ALEX2 and ALEX3, on the X chromosome. [provided by RefSeq, Jul 2008],similarity:Contains 3 ARM repeats.,tissue specificity:Expressed at high levels ovary, heart, testis, prostate, brain, spleen and colon. Expressed at very low levels in liver and thymus. Not expressed in peripheral blood leukocytes. Not or reduced expressed in lung, prostate, colon, pancreas and ovarian carcinomas.,

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of lysates from rat brain cells, using ARMCX1 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using ARMCX1 Polyclonal Antibody diluted at 1: 1000

产品名称: ARMCX1 Rabbit Polyclonal Antibody
产品货号: APRab07149



Western Blot analysis of mouse-brain cells using ARMCX1 Polyclonal Antibody diluted at 1: 1000

注意事项 (Note)

For research use only .