

产品名称: ARMCX3 Rabbit Polyclonal Antibody
产品货号: APRab07151



产品概述 (Summary)

产品名称 (Production Name)	ARMCX3 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ARMCX3
别名 (Alternative Names)	ARMCX3; ALEX3; BM-017; Armadillo repeat-containing X-linked protein 3; ARM protein lost in epithelial cancers on chromosome X 3; Protein ALEX3
基因 ID (Gene ID)	51566.0
蛋白 ID (SwissProt ID)	Q9UH62.The antiserum was produced against synthesized peptide derived from human ARMX3. AA range:291-340

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:5000-1:10000
蛋白分子量 (Molecular Weight)	38kDa

产品名称: ARMCX3 Rabbit Polyclonal Antibody
产品货号: APRab07151




研究背景 (Background)


armadillo repeat containing, X-linked 3(ARMCX3) Homo sapiens This gene encodes a member of the ALEX family of proteins which may play a role in tumor suppression. The encoded protein contains a potential N-terminal transmembrane domain and a single Armadillo (arm) repeat. Other proteins containing the arm repeat are involved in development, maintenance of tissue integrity, and tumorigenesis. This gene is closely localized with other family members on the X chromosome. Three transcript variants encoding the same protein have been identified for this gene. [provided by RefSeq, Jul 2008],similarity:Contains 3 ARM repeats.,

研究领域 (Research Area)

图片 (Image Data)




Immunofluorescence analysis of HepG2 cells, using ARMCX3 Antibody. The picture on the right is blocked with the synthesized peptide.




Immunohistochemistry analysis of paraffin-embedded human brain tissue, using ARMCX3 Antibody. The picture on the right is blocked with the synthesized peptide.

产品名称: ARMCX3 Rabbit Polyclonal Antibody
产品货号: APRab07151



Western blot analysis of lysates from K562 cells, using ARMCX3 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from COLO205 cells using ARMCX3 antibody.

注意事项 (Note)

For research use only .