

产品名称: ARP19 Rabbit Polyclonal Antibody
产品货号: APRab07157



产品概述 (Summary)

产品名称 (Production Name)	ARP19 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ARPP19
别名 (Alternative Names)	
基因 ID (Gene ID)	10776.0
蛋白 ID (SwissProt ID)	P56211.Synthesized peptide derived from part region of human protein AA range: 1-50

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	12kDa

研究背景 (Background)

The 19-kD cAMP-regulated phosphoprotein plays a role in regulating mitosis by inhibiting protein phosphatase-2A (PP2A;


产品名称: ARP19 Rabbit Polyclonal Antibody
产品货号: APRab07157



see MIM 176915) (summary by Gharbi-Ayachi et al., 2010 [PubMed 21164014]).[supplied by OMIM, Feb 2011],function:May enhance GAP-43 expression by binding to the NGF-regulatory region of its mRNA.,PTM:Phosphorylated by PKA.,similarity:Belongs to the endosulfine family.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from SW480 cells, primary antibody was diluted at 1:1000, 4°over night

注意事项 (Note)

For research use only .