

产品名称: Arrdc1 Rabbit Polyclonal Antibody
产品货号: APRab07165



产品概述 (Summary)

产品名称 (Production Name)	Arrdc1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, ELISA
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ARRDC1
别名 (Alternative Names)	ARRDC1; Arrestin domain-containing protein 1
基因 ID (Gene ID)	92714.0
蛋白 ID (SwissProt ID)	Q8N5I2. The antiserum was produced against synthesized peptide derived from human ARRD1. AA range: 167-216

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:100-1:300, ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	48kDa

研究背景 (Background)


产品名称: Arrdc1 Rabbit Polyclonal Antibody
产品货号: AP Rab07165




similarity:Belongs to the arrestin family.,similarity:Belongs to the arrestin family.,

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of lysates from Jurkat cells, using ARRDC1 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using Arrdc1 Polyclonal Antibody diluted at 1: 1000



Western Blot analysis of Jurkat cells using Arrdc1 Polyclonal Antibody diluted at 1: 1000

产品名称: Arrdc1 Rabbit Polyclonal Antibody
产品货号: APRab07165



注意事项 (Note)

For research use only .