

产品名称: ARRDC2 Rabbit Polyclonal Antibody
产品货号: APRab07166



产品概述 (Summary)

产品名称 (Production Name)	ARRDC2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ARRDC2
别名 (Alternative Names)	ARRDC2; PP2703; Arrestin domain-containing protein 2
基因 ID (Gene ID)	27106.0
蛋白 ID (SwissProt ID)	Q8TBH0.The antiserum was produced against synthesized peptide derived from human ARRDC2. AA range:141-190

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	44kDa

研究背景 (Background)


产品名称: ARRDC2 Rabbit Polyclonal Antibody
产品货号: AP Rab07166




similarity:Belongs to the arrestin family.,similarity:Belongs to the arrestin family.,

研究领域 (Research Area)


图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded human brain tissue, using ARRDC2 Antibody. The picture on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using ARRDC2 Polyclonal Antibody



Immunohistochemical analysis of paraffin-embedded Human brain. Antibody was diluted at 1:100 (4°,overnight) . High-pressure and temperature Tris-EDTA,pH8.0 was used for antigen retrieval. Negative contrl (right) obtained from antibody

产品名称: ARRDC2 Rabbit Polyclonal Antibody
产品货号: APRab07166



was pre-absorbed by immunogen peptide.

注意事项 (Note)

For research use only .