

产品名称: Arrdc4 Rabbit Polyclonal Antibody
产品货号: APRab07168



产品概述 (Summary)

产品名称 (Production Name)	Arrdc4 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ARRDC4
别名 (Alternative Names)	ARRDC4; Arrestin domain-containing protein 4
基因 ID (Gene ID)	91947.0
蛋白 ID (SwissProt ID)	Q8NCT1.The antiserum was produced against synthesized peptide derived from human ARRDC4. AA range:127-176

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	45kDa

研究背景 (Background)


产品名称: Arrdc4 Rabbit Polyclonal Antibody
产品货号: AP Rab07168



similarity:Belongs to the arrestin family.,similarity:Belongs to the arrestin family.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from COLO cells, using ARRDC4 Antibody. The lane on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .