

产品名称: Atg16L2 Rabbit Polyclonal Antibody
产品货号: APRab07287



产品概述 (Summary)

产品名称 (Production Name)	Atg16L2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ATG16L2
别名 (Alternative Names)	ATG16L2; WDR80; Autophagy-related protein 16-2; APG16-like 2; WD repeat-containing protein 80
基因 ID (Gene ID)	89849.0
蛋白 ID (SwissProt ID)	Q8NAA4.The antiserum was produced against synthesized peptide derived from human ATG16L2. AA range:241-290

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	70kDa

产品名称: Atg16L2 Rabbit Polyclonal Antibody
产品货号: APRab07287




研究背景 (Background)


function: May play a role in autophagy., similarity: Belongs to the WD repeat ATG16 family., similarity: Contains 7 WD repeats., function: May play a role in autophagy., similarity: Belongs to the WD repeat ATG16 family., similarity: Contains 7 WD repeats.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of ATG16L2 Antibody. The lane on the right is blocked with the ATG16L2 peptide.



Western blot analysis of the lysates from RAW264.7 cells using ATG16L2 antibody.

注意事项 (Note)

For research use only .