

产品名称: Atg9a Rabbit Polyclonal Antibody
产品货号: APRab07302



产品概述 (Summary)

产品名称 (Production Name)	Atg9a Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ATG9A APG9L1
别名 (Alternative Names)	ATG9 autophagy related 9 homolog A (S. cerevisiae)
基因 ID (Gene ID)	79065.0
蛋白 ID (SwissProt ID)	Q7Z3C6.Synthetic Peptide of Atg9a

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	85kDa

研究背景 (Background)

function:Plays a role in autophagy.,similarity:Belongs to the ATG9 family.,function:Plays a role in


产品名称: Atg9a Rabbit Polyclonal Antibody
产品货号: APRab07302



autophagy.,similarity:Belongs to the ATG9 family.,

研究领域 (Research Area)

图片 (Image Data)



Western Blot analysis of 293T 3T3 cells using Atg9a Polyclonal Antibody diluted at 1:800. Secondary antibody was diluted at 1:20000

注意事项 (Note)

For research use only .