

产品名称: ATP5S Rabbit Polyclonal Antibody
产品货号: APRab07340



产品概述 (Summary)

产品名称 (Production Name)	ATP5S Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ATP5S
别名 (Alternative Names)	ATP5S; ATPW; ATP synthase subunit s; mitochondrial; ATP synthase-coupling factor B; FB; Mitochondrial ATP synthase regulatory component factor B
基因 ID (Gene ID)	27109.0
蛋白 ID (SwissProt ID)	Q99766.The antiserum was produced against synthesized peptide derived from human ATP5S. AA range:21-70

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	23kDa

产品名称: ATP5S Rabbit Polyclonal Antibody
产品货号: APRab07340





研究背景 (Background)


This gene encodes a subunit of mitochondrial ATP synthase. Mitochondrial ATP synthase catalyzes ATP synthesis, utilizing an electrochemical gradient of protons across the inner membrane during oxidative phosphorylation. ATP synthase is composed of two linked multi-subunit complexes: the soluble catalytic core, F1, and the membrane-spanning component, Fo, comprising the proton channel. This gene encodes the subunit s, also known as factor B, of the proton channel. This subunit is necessary for the energy transduction activity of the ATP synthase complexes. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2008],caution:It is uncertain whether Met-1 or Met-16 is the initiator.,function:Involved in regulation of mitochondrial membrane ATP synthase. Necessary for H(+) conduction of ATP synthase.,similarity:Belongs to the ATP synthase subunit s family.,subunit:Monomer. Associates with ATP synthase.,

研究领域 (Research Area)

图片 (Image Data)



产品名称: ATP5S Rabbit Polyclonal Antibody
产品货号: APRab07340



Western Blot analysis of various cells using ATP5S Polyclonal Antibody

注意事项 (Note)

For research use only .