

产品名称: Aven Rabbit Polyclonal Antibody
产品货号: APRab07377



产品概述 (Summary)

产品名称 (Production Name)	Aven Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	AVEN
别名 (Alternative Names)	AVEN; Cell death regulator Aven
基因 ID (Gene ID)	57099.0
蛋白 ID (SwissProt ID)	Q9NQS1.The antiserum was produced against synthesized peptide derived from human AVEN. AA range:301-350

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	55kDa

研究背景 (Background)

产品名称: Aven Rabbit Polyclonal Antibody
产品货号: AP Rab07377





function: Protects against apoptosis mediated by Apaf-1, subcellular location: Associated with intracellular membranes, subunit: Binds Apaf-1, BCL-2 and BAD (Bcl-xl), tissue specificity: Highly expressed in testis, ovary, thymus, prostate, spleen, small intestine, colon, heart, skeletal muscle, liver, kidney and pancreas, function: Protects against apoptosis mediated by Apaf-1, subcellular location: Associated with intracellular membranes, subunit: Binds Apaf-1, BCL-2 and BAD (Bcl-xl), tissue specificity: Highly expressed in testis, ovary, thymus, prostate, spleen, small intestine, colon, heart, skeletal muscle, liver, kidney and pancreas.

研究领域 (Research Area)

Cell Biology

图片 (Image Data)



注意事项 (Note)

For research use only .