

产品名称: BATF3 Rabbit Polyclonal Antibody
产品货号: APRab07472



产品概述 (Summary)

产品名称 (Production Name)	BATF3 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	BATF3 BATF3; SNFT; Basic leucine zipper transcriptional factor ATF-like 3; B-ATF-3; 21
别名 (Alternative Names)	kDa small nuclear factor isolated from T-cells; Jun dimerization protein p21SNFT
基因 ID (Gene ID)	55509.0
蛋白 ID (SwissProt ID)	Q9NR55.The antiserum was produced against synthesized peptide derived from the Internal region of human BATF3. AA range:1-50

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	15kDa

产品名称: BATF3 Rabbit Polyclonal Antibody
产品货号: APRab07472




研究背景 (Background)


This gene encodes a member of the basic leucine zipper protein family. The encoded protein functions as a transcriptional repressor when heterodimerizing with JUN. The protein may play a role in repression of interleukin-2 and matrix metalloproteinase-1 transcription.[provided by RefSeq, Feb 2009],function:Functions as a negative regulator of AP-1-mediated transcription by heterodimerizing with JUN and binding DNA at 12-O-tetradecanoylphorbol-13-acetate response elements (TRE) (consensus: 5'-TGA[CG]TCA-3'). Represses IL2 and MMP1 promoter activities.,similarity:Belongs to the bZIP family.,similarity:Contains 1 bZIP domain.,subunit:Interacts with JUN.,

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of mouse-lung mouse-spleen lysis using BATF3 antibody. Antibody was diluted at 1:2000. Secondary antibody was diluted at 1:20000 cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Invent biotech, MN, USA) .



Immunohistochemical analysis of paraffin-embedded human-breast-cancer, antibody was diluted at 1:200

产品名称: BATF3 Rabbit Polyclonal Antibody
产品货号: APRab07472



Immunohistochemical analysis of paraffin-embedded human-breast-cancer, antibody was diluted at 1:200

注意事项 (Note)

For research use only .