

产品名称: BD1L1 Rabbit Polyclonal Antibody
产品货号: APRab07522



产品概述 (Summary)

产品名称 (Production Name)	BD1L1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC/IF
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	BOD1L1
别名 (Alternative Names)	BOD1L FAM44A KIAA1327
基因 ID (Gene ID)	259282.0
蛋白 ID (SwissProt ID)	Q8NFC6.Synthesized peptide derived from human protein . at AA range: 1280-1360

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:50-1:300, ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	335kDa

研究背景 (Background)

Component of the fork protection machinery required to protect stalled/damaged replication forks from uncontrolled



研究领域 (Research Area)

注意事项 (Note)

For research use only .