

产品名称: BEGAIN Rabbit Polyclonal Antibody
产品货号: APRab07532



产品概述 (Summary)

产品名称 (Production Name)	BEGAIN Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	BEGAIN
别名 (Alternative Names)	BEGAIN; KIAA1446; Brain-enriched guanylate kinase-associated protein
基因 ID (Gene ID)	57596.0
蛋白 ID (SwissProt ID)	Q9BUH8.The antiserum was produced against synthesized peptide derived from human BEGIN. AA range:511-560

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:5000-1:10000
蛋白分子量 (Molecular Weight)	65kDa

研究背景 (Background)


产品名称: BEGAIN Rabbit Polyclonal Antibody
产品货号: APRab07532



function:May sustain the structure of the postsynaptic density (PSD).,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,subunit:Interacts with DLG4 and DLGAP1 and forms a ternary complex.,function:May sustain the structure of the postsynaptic density (PSD).,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,subunit:Interacts with DLG4 and DLGAP1 and forms a ternary complex.,

研究领域 (Research Area)

图片 (Image Data)



Immunofluorescence analysis of HeLa cells, using BEGAIN Antibody. The picture on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .