

产品名称: BNIP-2 Rabbit Polyclonal Antibody
产品货号: APRab07616



产品概述 (Summary)

产品名称 (Production Name)	BNIP-2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	BNIP2
别名 (Alternative Names)	BNIP2; NIP2; BCL2/adenovirus E1B 19 kDa protein-interacting protein 2
基因 ID (Gene ID)	663.0
蛋白 ID (SwissProt ID)	Q12982.The antiserum was produced against synthesized peptide derived from human BNIP-2. AA range:73-122

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	36kDa

研究背景 (Background)


产品名称: BNIP-2 Rabbit Polyclonal Antibody
产品货号: APRab07616




This gene is a member of the BCL2/adenovirus E1B 19 kd-interacting protein (BNIP) family. It interacts with the E1B 19 kDa protein, which protects cells from virally-induced cell death. The encoded protein also interacts with E1B 19 kDa-like sequences of BCL2, another apoptotic protector. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Mar 2016],function:Implicated in the suppression of cell death. Interacts with the BCL-2 and adenovirus E1B 19 kDa proteins.,similarity:Contains 1 CRAL-TRIO domain.,subcellular location:Localizes to the nuclear envelope region and to other cytoplasmic structures.,

研究领域 (Research Area)

图片 (Image Data)




Immunohistochemistry analysis of BNIP-2 antibody in paraffin-embedded human breast carcinoma tissue.




Western blot analysis of lysate from Jurkat cells, using BNIP-2 antibody.

产品名称: BNIP-2 Rabbit Polyclonal Antibody
产品货号: APRab07616



Western Blot analysis of various cells using BNIP-2 Polyclonal Antibody



Western Blot analysis of K562 cells using BNIP-2 Polyclonal Antibody

注意事项 (Note)

For research use only .