

产品名称: BNIP-2 Rabbit Polyclonal Antibody
产品货号: APRab07617



产品概述 (Summary)

产品名称 (Production Name)	BNIP-2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, ELISA, FC
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	BNIP2
别名 (Alternative Names)	BNIP2; NIP2; BCL2/adenovirus E1B 19 kDa protein-interacting protein 2
基因 ID (Gene ID)	663.0
蛋白 ID (SwissProt ID)	Q12982. The antiserum was produced against synthesized peptide derived from human BNIP2. AA range: 265-314

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:200-1:1000, ELISA 1:20000-1:40000, FC 1:200-1:400
蛋白分子量 (Molecular Weight)	36kDa

研究背景 (Background)


产品名称: BNIP-2 Rabbit Polyclonal Antibody
产品货号: AP Rab07617



This gene is a member of the BCL2/adenovirus E1B 19 kd-interacting protein (BNIP) family. It interacts with the E1B 19 kDa protein, which protects cells from virally-induced cell death. The encoded protein also interacts with E1B 19 kDa-like sequences of BCL2, another apoptotic protector. Alternate splicing results in multiple transcript variants. [provided by RefSeq, Mar 2016],function:Implicated in the suppression of cell death. Interacts with the BCL-2 and adenovirus E1B 19 kDa proteins.,similarity:Contains 1 CRAL-TRIO domain.,subcellular location:Localizes to the nuclear envelope region and to other cytoplasmic structures.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from Jurkat cells, using BNIP2 Antibody. The lane on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .