

产品名称: BRMS-1 Rabbit Polyclonal Antibody
产品货号: APRab07660



产品概述 (Summary)

产品名称 (Production Name)	BRMS-1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	BRMS1
别名 (Alternative Names)	BRMS1; Breast cancer metastasis-suppressor 1
基因 ID (Gene ID)	25855.0
蛋白 ID (SwissProt ID)	Q9HCU9.The antiserum was produced against synthesized peptide derived from the C-terminal region of human BRMS1. AA range:197-246

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	28kDa

研究背景 (Background)


产品名称: BRMS-1 Rabbit Polyclonal Antibody
产品货号: APRab07660




This gene reduces the metastatic potential, but not the tumorigenicity, of human breast cancer and melanoma cell lines. The protein encoded by this gene localizes primarily to the nucleus and is a component of the mSin3a family of histone deacetylase complexes (HDAC). The protein contains two coiled-coil motifs and several imperfect leucine zipper motifs. Alternative splicing results in two transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008],function:May be a mediator of metastasis suppression in breast carcinoma.,similarity:Belongs to the BRMS1 family.,tissue specificity:Expression levels are higher in term placentas than in early placentas. Low levels of expression observed in normal pregnancies and in molar pregnancies.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysate from HT29 cells, using BRMS1 Antibody.



Western Blot analysis of HT29, K562, BT474 cells using BRMS-1 Polyclonal Antibody.. Secondary antibody was diluted at 1:20000

注意事项 (Note)

For research use only .