

产品名称: Brp16 Rabbit Polyclonal Antibody
产品货号: AP Rab07663



产品概述 (Summary)

产品名称 (Production Name)	Brp16 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FAM203A
别名 (Alternative Names)	FAM203A; BRP16; C8orf30A; Protein FAM203A; Brain protein 16
基因 ID (Gene ID)	51236.0
蛋白 ID (SwissProt ID)	Q9BTY7.The antiserum was produced against synthesized peptide derived from human BRP16. AA range:251-300

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	40kDa

研究背景 (Background)


产品名称: Brp16 Rabbit Polyclonal Antibody
产品货号: APRab07663




similarity:Belongs to the UPF0570 (BRP16/HGH1) family.,similarity:Belongs to the UPF0570 (BRP16/HGH1) family.,

研究领域 (Research Area)

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue, using BRP16 Antibody. The picture on the right is blocked with the synthesized peptide.



Western blot analysis of lysates from COLO cells, using BRP16 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using Brp16 Polyclonal Antibody

产品名称: Brp16 Rabbit Polyclonal Antibody
产品货号: APRab07663



注意事项 (Note)

For research use only .