

产品名称: BS69 Rabbit Polyclonal Antibody
产品货号: APRab07672



产品概述 (Summary)

产品名称 (Production Name)	BS69 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ZMYND11
别名 (Alternative Names)	ZMYND11; BS69; Zinc finger MYND domain-containing protein 11; Adenovirus 5 E1A-binding protein; Protein BS69
基因 ID (Gene ID)	10771.0
蛋白 ID (SwissProt ID)	Q15326. The antiserum was produced against synthesized peptide derived from human ZMY11. AA range: 111-160

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:300, ICC/IF 1:50-1:200, ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	

产品名称: BS69 Rabbit Polyclonal Antibody
产品货号: APRab07672




研究背景 (Background)

The protein encoded by this gene was first identified by its ability to bind the adenovirus E1A protein. The protein localizes to the nucleus. It functions as a transcriptional repressor, and expression of E1A inhibits this repression. Alternatively spliced transcript variants encoding different isoforms have been identified. [provided by RefSeq, Jul 2008],function: Binds to the transactivation domain of the adenovirus type 5 E1A 32 kDa protein (289R) and inhibits its transactivating activity. May act as tumor suppressor through suppression of adenovirus replication.,PTM: Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity: Contains 1 bromo domain.,similarity: Contains 1 MYND-type zinc finger.,similarity: Contains 1 PHD-type zinc finger.,similarity: Contains 1 PWWP domain.,

研究领域 (Research Area)

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using ZMY11 Antibody. The picture on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .