

产品名称: BTF Rabbit Polyclonal Antibody
产品货号: APRab07685



产品概述 (Summary)

产品名称 (Production Name)	BTF Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	BCLAF1
别名 (Alternative Names)	BCLAF1; BTF; KIAA0164; Bcl-2-associated transcription factor 1; Btf
基因 ID (Gene ID)	9774.0
蛋白 ID (SwissProt ID)	Q9NYF8.The antiserum was produced against synthesized peptide derived from human BCLAF1. AA range:281-330

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	106kDa

研究背景 (Background)


产品名称: BTF Rabbit Polyclonal Antibody
产品货号: APRab07685




This gene encodes a transcriptional repressor that interacts with several members of the BCL2 family of proteins. Overexpression of this protein induces apoptosis, which can be suppressed by co-expression of BCL2 proteins. The protein localizes to dot-like structures throughout the nucleus, and redistributes to a zone near the nuclear envelope in cells undergoing apoptosis. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008],function:Death-promoting transcriptional repressor.,sequence caution:Contaminating sequence. Potential poly-A sequence.,subunit:Interacts with Bcl-2 related proteins, EMD, with the adenovirus E1B 19 kDa protein and with DNA.,tissue specificity:Ubiquitous.,

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of the lysates from HeLa cells using BCLAF1 antibody.



Western Blot analysis of various cells using BTF Polyclonal Antibody

产品名称: BTF Rabbit Polyclonal Antibody
产品货号: APRab07685



Western Blot analysis of A549 cells using BTF Polyclonal Antibody

注意事项 (Note)

For research use only .