

产品名称: Bub3 Rabbit Polyclonal Antibody
产品货号: APRab07695



产品概述 (Summary)

产品名称 (Production Name)	Bub3 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	BUB3
别名 (Alternative Names)	BUB3; Mitotic checkpoint protein BUB3
基因 ID (Gene ID)	9184.0
蛋白 ID (SwissProt ID)	O43684.The antiserum was produced against synthesized peptide derived from human BUB3. AA range:279-328

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:5000-1:10000
蛋白分子量 (Molecular Weight)	40kDa

研究背景 (Background)

产品名称: Bub3 Rabbit Polyclonal Antibody
产品货号: AP Rab07695




This gene encodes a protein involved in spindle checkpoint function. The encoded protein contains four WD repeat domains and has sequence similarity with the yeast BUB3 protein. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. [provided by RefSeq, Jul 2008],function:Required for kinetochore localization of BUB1.,similarity:Belongs to the WD repeat BUB3 family.,similarity:Contains 5 WD repeats.,subunit:Interacts with BUB1 and BUBR1.,


研究领域 (Research Area)

Cell_Cycle_G1S;Cell_Cycle_G2M_DNA;

图片 (Image Data)




Western blot analysis of lysates from HeLa cells, Jurkat cells, and 293 cells, using BUB3 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using Bub3 Polyclonal Antibody

产品名称: Bub3 Rabbit Polyclonal Antibody
产品货号: AP Rab07695



Western Blot analysis of 293 cells using Bub3 Polyclonal Antibody

注意事项 (Note)

For research use only .