

产品名称: Cadherin-24 Rabbit Polyclonal Antibody
产品货号: AP Rab07833



产品概述 (Summary)

产品名称 (Production Name)	Cadherin-24 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, ELISA
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CDH24
别名 (Alternative Names)	CDH24; CDH11L; Cadherin-24
基因 ID (Gene ID)	64403.0
蛋白 ID (SwissProt ID)	Q86UP0. The antiserum was produced against synthesized peptide derived from human CDH24. AA range: 81-130

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:200-1:1000, ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	

研究背景 (Background)

产品名称: Cadherin-24 Rabbit Polyclonal Antibody
产品货号: APRab07833




function: Cadherins are calcium dependent cell adhesion proteins. They preferentially interact with themselves in a homophilic manner in connecting cells; cadherins may thus contribute to the sorting of heterogeneous cell types. Cadherin-24 mediate strong cell-cell adhesion. similarity: Contains 5 cadherin domains. subunit: Associates with alpha-, beta- and delta-catenins. function: Cadherins are calcium dependent cell adhesion proteins. They preferentially interact with themselves in a homophilic manner in connecting cells; cadherins may thus contribute to the sorting of heterogeneous cell types. Cadherin-24 mediate strong cell-cell adhesion. similarity: Contains 5 cadherin domains. subunit: Associates with alpha-, beta- and delta-catenins.

研究领域 (Research Area)

Adherens_Junction

图片 (Image Data)



Immunofluorescence analysis of HepG2 cells, using CDH24 Antibody. The picture on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .