

产品名称: Cadherin-26 Rabbit Polyclonal Antibody
产品货号: APRab07834



产品概述 (Summary)

产品名称 (Production Name)	Cadherin-26 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CDH26
别名 (Alternative Names)	CDH26; Cadherin-like protein 26; Cadherin-like protein VR20
基因 ID (Gene ID)	60437.0
蛋白 ID (SwissProt ID)	Q8IXH8. The antiserum was produced against synthesized peptide derived from human CDH26. AA range: 11-60

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:200-1:1000, ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	95kDa

研究背景 (Background)

产品名称: Cadherin-26 Rabbit Polyclonal Antibody
产品货号: APRab07834




Cadherins are a family of adhesion molecules that mediate Ca^{2+} -dependent cell-cell adhesion in all solid tissues and modulate a wide variety of processes, including cell polarization and migration. Cadherin domains occur as repeats in the extracellular region and are thought to contribute to the sorting of heterogeneous cell types and the maintenance of orderly structures such as epithelium. This gene encodes a cadherin domain-containing protein whose specific function has not yet been determined. Alternative splicing occurs at this locus and two transcript variants, encoding distinct proteins, have been identified. [provided by RefSeq, Jul 2008],function:Cadherins are calcium dependent cell adhesion proteins. They preferentially interact with themselves in a homophilic manner in connecting cells; cadherins may thus contribute to the sorting of heterogeneous cell types.,similarity:Contains 4 cadherin domains.,

研究领域 (Research Area)

Adherens_Junction

图片 (Image Data)



Immunofluorescence analysis of MCF7 cells, using CDH26 Antibody. The picture on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using Cadherin-26 Polyclonal Antibody

注意事项 (Note)

For research use only .

产品名称: Cadherin-26 Rabbit Polyclonal Antibody
产品货号: APRab07834

