

产品名称: CAGE-1 Rabbit Polyclonal Antibody
产品货号: APRab07840



产品概述 (Summary)

产品名称 (Production Name)	CAGE-1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CAGE1
别名 (Alternative Names)	CAGE1; CTAG3; Cancer-associated gene 1 protein; CAGE-1; Cancer/testis antigen 3; CT3
基因 ID (Gene ID)	285782.0
蛋白 ID (SwissProt ID)	Q8TC20.The antiserum was produced against synthesized peptide derived from human CAGE1. AA range:711-760

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	90kDa

产品名称: CAGE-1 Rabbit Polyclonal Antibody
产品货号: APRab07840




研究背景 (Background)


tissue specificity: Testis-specific expression in normal tissues, but wide expression among cancer tissues and cell lines., tissue specificity: Testis-specific expression in normal tissues, but wide expression among cancer tissues and cell lines.,

研究领域 (Research Area)

图片 (Image Data)




Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using CAGE1 Antibody. The picture on the right is blocked with the synthesized peptide.



Western blot analysis of lysates from HT-29 and Jurkat cells, using CAGE1 Antibody. The lane on the right is blocked with the synthesized peptide.

产品名称: CAGE-1 Rabbit Polyclonal Antibody
产品货号: APRab07840



Western blot analysis of the lysates from HepG2 cells using CAGE1 antibody.



Western Blot analysis of various cells using CAGE-1 Polyclonal Antibody diluted at 1: 2000

注意事项 (Note)

For research use only .