

产品名称: CARD 6 Rabbit Polyclonal Antibody
产品货号: APRab07926



产品概述 (Summary)

产品名称 (Production Name)	CARD 6 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CARD6
别名 (Alternative Names)	CARD6; Caspase recruitment domain-containing protein 6
基因 ID (Gene ID)	84674.0
蛋白 ID (SwissProt ID)	Q9BX69.The antiserum was produced against synthesized peptide derived from human CARD6. AA range:1-50

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	116kDa

研究背景 (Background)

产品名称: CARD 6 Rabbit Polyclonal Antibody
产品货号: AP Rab07926




This gene encodes a protein that contains a caspase recruitment domain (CARD), an antiparallel six-helical bundle that mediates homotypic protein-protein interactions. The encoded protein is a microtubule-associated protein that has been shown to interact with receptor-interacting protein kinases and positively modulate signal transduction pathways converging on activation of the inducible transcription factor NF- κ B. [provided by RefSeq, Jul 2008],function:May be involved in apoptosis.,similarity:Contains 1 CARD domain.,


研究领域 (Research Area)

NOD-like receptor;

图片 (Image Data)




Immunofluorescence analysis of HeLa cells, using CARD6 Antibody. The picture on the right is blocked with the synthesized peptide.




Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using CARD6 Antibody. The picture on the right is blocked with the synthesized peptide.

产品名称: CARD 6 Rabbit Polyclonal Antibody
产品货号: APRab07926



Western blot analysis of lysates from HeLa cells, using CARD6 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using CARD 6 Polyclonal Antibody

注意事项 (Note)

For research use only .