

产品名称: CCDC4 Rabbit Polyclonal Antibody
产品货号: APRab08119



产品概述 (Summary)

产品名称 (Production Name)	CCDC4 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	BEND4
别名 (Alternative Names)	BEND4; CCDC4; BEN domain-containing protein 4; Coiled-coil domain-containing protein 4
基因 ID (Gene ID)	389206.0
蛋白 ID (SwissProt ID)	Q6ZU67.The antiserum was produced against synthesized peptide derived from human BEND4. AA range:371-420

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	57kDa

产品名称: CCDC4 Rabbit Polyclonal Antibody
产品货号: APRab08119




研究背景 (Background)

PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 1 BEN domain.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 1 BEN domain.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of BEND4 Antibody. The lane on the right is blocked with the BEND4 peptide.



Western Blot analysis of various cells using CCDC4 Polyclonal Antibody

注意事项 (Note)

For research use only .