

产品名称: CCDC45 Rabbit Polyclonal Antibody
产品货号: APRab08120



产品概述 (Summary)

产品名称 (Production Name)	CCDC45 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CEP95
别名 (Alternative Names)	CEP95; CCDC45; CEP45; Centrosomal protein of 95 kDa; Coiled-coil domain-containing protein 45
基因 ID (Gene ID)	90799.0
蛋白 ID (SwissProt ID)	Q96GE4.The antiserum was produced against synthesized peptide derived from human CCDC45. AA range:591-640

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	95kDa

产品名称: CCDC45 Rabbit Polyclonal Antibody
产品货号: APRab08120




研究背景 (Background)


PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of the lysates from Jurkat cells using CCDC45 antibody.



Western blot analysis of various cells using CCDC45 Polyclonal Antibody



产品名称: CCDC45 Rabbit Polyclonal Antibody
产品货号: APRab08120



注意事项 (Note)

For research use only .