

产品名称: CCDC99 Rabbit Polyclonal Antibody
产品货号: APRab08124



产品概述 (Summary)

产品名称 (Production Name)	CCDC99 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CCDC99
别名 (Alternative Names)	CCDC99; Protein Spindly; hSpindly; Arsenite-related gene 1 protein; Coiled-coil domain-containing protein 99; Rhabdomyosarcoma antigen MU-RMS-40.4A
基因 ID (Gene ID)	54908.0
蛋白 ID (SwissProt ID)	Q96EA4.The antiserum was produced against synthesized peptide derived from human CCDC99. AA range:541-590

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	70kDa

产品名称: CCDC99 Rabbit Polyclonal Antibody
产品货号: APRab08124




研究背景 (Background)

This gene encodes a coiled-coil domain-containing protein that functions in mitotic spindle formation and chromosome segregation. The encoded protein plays a role in coordinating microtubule attachment by promoting recruitment of dynein proteins, and in mitotic checkpoint signaling. [provided by RefSeq, Jul 2016], mitotic sister chromatid segregation, M phase of mitotic cell cycle, establishment of mitotic spindle orientation, microtubule cytoskeleton organization, mitotic cell cycle, M phase, nuclear division, sister chromatid segregation, cell morphogenesis, cytoskeleton organization, microtubule-based process, cell cycle, chromosome segregation, mitosis, mitotic metaphase plate congression, establishment or maintenance of cell polarity, protein localization, cell cycle process, cell cycle phase, establishment of cell polarity, cellular component morphogenesis, protein localization to kinetochore, cellular protein localization, establishment of mitotic spindle localization, organelle fission, chromosome localization, chromosome organization, establishment of spindle localization, establishment of spindle orientation, establishment of chromosome localization, metaphase plate congression, organelle localization, spindle localization, establishment of organelle localization, cellular macromolecule localization,


研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of CCDC99 Antibody. The lane on the right is blocked with the CCDC99 peptide.

产品名称: CCDC99 Rabbit Polyclonal Antibody
产品货号: APRab08124



Western blot analysis of the lysates from K562 cells using CCDC99 antibody.



Western Blot analysis of 3T3 cells using CCDC99 Polyclonal Antibody diluted at 1: 1000

注意事项 (Note)

For research use only .