

产品名称: CCNI Rabbit Polyclonal Antibody
产品货号: APRab08149



产品概述 (Summary)

产品名称 (Production Name)	CCNI Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CCNI
别名 (Alternative Names)	
基因 ID (Gene ID)	10983.0
蛋白 ID (SwissProt ID)	Q14094.Synthesized peptide derived from human protein . at AA range: 30-110

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	41kDa

研究背景 (Background)

The protein encoded by this gene belongs to the highly conserved cyclin family, whose members are characterized by a


产品名称: CCNI Rabbit Polyclonal Antibody
产品货号: APRab08149



dramatic periodicity in protein abundance through the cell cycle. Cyclins function as regulators of CDK kinases. Different cyclins exhibit distinct expression and degradation patterns which contribute to the temporal coordination of each mitotic event. This cyclin shows the highest similarity with cyclin G. The transcript of this gene was found to be expressed constantly during cell cycle progression. The function of this cyclin has not yet been determined. [provided by RefSeq, Jul 2008],developmental stage:Expression is independent of the cell cycle in lung fibroblasts.,similarity:Belongs to the cyclin family.,tissue specificity:Highest levels in adult heart, brain and skeletal muscle. Lower levels in adult placenta, lung, kidney and pancreas. Also high levels in fetal brain and lower levels in fetal lung, liver and kidney. Also abundant in testis and thyroid.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from SW480 cells, primary antibody was diluted at 1:1000, 4°over night

注意事项 (Note)

For research use only .