

**产品名称: CCNK Rabbit Polyclonal Antibody**  
**产品货号: APRab08150**



## 产品概述 (Summary)

|                        |                                 |
|------------------------|---------------------------------|
| 产品名称 (Production Name) | CCNK Rabbit Polyclonal Antibody |
| 描述 (Description)       | Rabbit polyclonal Antibody      |
| 宿主 (Host)              | Rabbit                          |
| 应用 (Application)       | WB,ELISA                        |
| 种属反应性 (Reactivity)     | Human,Mouse                     |

## 产品性能 (Performance)

|                     |                                                                                          |
|---------------------|------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                             |
| 修饰 (Modification)   | Unmodified                                                                               |
| 同种型 (Isotype)       | IgG                                                                                      |
| 克隆 (Clonality)      | Polyclonal                                                                               |
| 形式 (Form)           | Liquid                                                                                   |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles. |
| 储存溶液 (Buffer)       | Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.                |
| 纯化方式 (Purification) | Affinity purification                                                                    |

## 免疫原信息 (Immunogen)

|                        |                                                                             |
|------------------------|-----------------------------------------------------------------------------|
| 基因名 (Gene Name)        | CCNK CPR4                                                                   |
| 别名 (Alternative Names) |                                                                             |
| 基因 ID (Gene ID)        | 8812.0                                                                      |
| 蛋白 ID (SwissProt ID)   | O75909.Synthesized peptide derived from human protein . at AA range: 80-160 |

## 产品应用 (Application)

|                          |                                      |
|--------------------------|--------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:2000,ELISA 1:5000-1:20000 |
| 蛋白分子量 (Molecular Weight) | 63kDa                                |

## 研究背景 (Background)

The protein encoded by this gene is a member of the transcription cyclin family. These cyclins may regulate transcription


**产品名称: CCNK Rabbit Polyclonal Antibody**  
**产品货号: APRab08150**



through their association with and activation of cyclin-dependent kinases (CDK) that phosphorylate the C-terminal domain (CTD) of the large subunit of RNA polymerase II. This gene product may play a dual role in regulating CDK and RNA polymerase II activities. [provided by RefSeq, Jul 2008],function:May play a role in transcriptional regulation. In vitro, is associated with a kinase activity toward both RNA polymerase II C-terminal domain and CDK2 (CAK).,similarity:Belongs to the cyclin family. Cyclin C subfamily.,subunit:Part of a cyclin-kinase pair in the RNA polymerase II holoenzyme. Binds to CDK9.,tissue specificity:Ubiquitously expressed. Highest levels in testis.,

## 研究领域 (Research Area)

## 图片 (Image Data)



Western Blot analysis of HEK293 lysis, using primary antibody at 1:1000 dilution. Secondary antibody was diluted at 1:10000

## 注意事项 (Note)

For research use only .