

产品名称: CCP2 Rabbit Polyclonal Antibody
产品货号: APRab08155



产品概述 (Summary)

产品名称 (Production Name)	CCP2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	AGBL2
别名 (Alternative Names)	AGBL2; CCP2; Cytosolic carboxypeptidase 2; ATP/GTP-binding protein-like 2
基因 ID (Gene ID)	79841.0
蛋白 ID (SwissProt ID)	Q5U5Z8.The antiserum was produced against synthesized peptide derived from human CBCP2. AA range:731-780

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	104kDa

研究背景 (Background)


产品名称: CCP2 Rabbit Polyclonal Antibody
产品货号: APRab08155




cofactor: Binds 1 zinc ion per subunit., function: May play a role in the processing of tubulin., sequence caution: Translated as Lys., similarity: Belongs to the peptidase M14 family., cofactor: Binds 1 zinc ion per subunit., function: May play a role in the processing of tubulin., sequence caution: Translated as Lys., similarity: Belongs to the peptidase M14 family.,

研究领域 (Research Area)

图片 (Image Data)




Immunohistochemistry analysis of paraffin-embedded human heart tissue, using CBCP2 Antibody. The picture on the right is blocked with the synthesized peptide.




Western blot analysis of lysates from HUVEC cells, using CBCP2 Antibody. The lane on the right is blocked with the synthesized peptide.

产品名称: CCP2 Rabbit Polyclonal Antibody
产品货号: APRab08155



Western blot analysis of the lysates from Jurkat cells using CBCP2 antibody.



Western Blot analysis of various cells using CCP2 Polyclonal Antibody diluted at 1: 2000

注意事项 (Note)

For research use only .