

产品名称: CCRK Rabbit Polyclonal Antibody
产品货号: APRab08165



产品概述 (Summary)

产品名称 (Production Name)	CCRK Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CDK20 CDK20; CCRK; CDCH; Cyclin-dependent kinase 20; CDK-activating kinase p42;
别名 (Alternative Names)	CAK-kinase p42; Cell cycle-related kinase; Cell division protein kinase 20; Cyclin-dependent protein kinase H; Cyclin-kinase-activating kinase p42
基因 ID (Gene ID)	23552.0
蛋白 ID (SwissProt ID)	Q8IZL9.The antiserum was produced against synthesized peptide derived from human CCRK. AA range:31-80

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	39kDa

产品名称: CCRK Rabbit Polyclonal Antibody
产品货号: APRab08165




研究背景 (Background)


The protein encoded by this gene contains a kinase domain most closely related to the cyclin-dependent protein kinases. The encoded kinase may activate cyclin-dependent kinase 2 and is involved in cell growth. Alternatively spliced transcript variants encoding distinct isoforms have been reported. [provided by RefSeq, Dec 2009],similarity:Belongs to the protein kinase superfamily.,

研究领域 (Research Area)

图片 (Image Data)




Immunofluorescence analysis of HUVEC cells, using CCRK Antibody. The picture on the right is blocked with the synthesized peptide.



Western blot analysis of lysates from RAW264.7 cells, using CCRK Antibody. The lane on the right is blocked with the synthesized peptide.

产品名称: CCRK Rabbit Polyclonal Antibody
产品货号: APRab08165



Western blot analysis of the lysates from RAW264.7 cells using CCRK antibody.

注意事项 (Note)

For research use only .