

产品名称: CD300c Rabbit Polyclonal Antibody
产品货号: APRab08341



产品概述 (Summary)

产品名称 (Production Name)	CD300c Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD300C CD300C; CMRF35; CMRF35A; CMRF35A1; IGSF16; CMRF35-like molecule 6;
别名 (Alternative Names)	CLM-6; CD300 antigen-like family member C; CMRF35-A1; CMRF-35; Immunoglobulin superfamily member 16; IgSF16; CD300c
基因 ID (Gene ID)	10871.0
蛋白 ID (SwissProt ID)	Q08708.The antiserum was produced against synthesized peptide derived from the Internal region of human CD300C. AA range:131-180

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	37kDa

产品名称: CD300c Rabbit Polyclonal Antibody
产品货号: APRab08341




研究背景 (Background)


CD300c molecule(CD300C) Homo sapiens The CMRF35 antigen, which was identified by reactivity with a monoclonal antibody, is present on monocytes, neutrophils, and some T and B lymphocytes (Jackson et al., 1992 [PubMed 1349532]). [supplied by OMIM, Mar 2008],similarity:Belongs to the CD300 family.,similarity:Contains 1 Ig-like V-type (immunoglobulin-like) domain.,tissue specificity:Present on the surface of monocytes, neutrophils, a proportion of peripheral blood T- and B-lymphocytes and lymphocytic cell lines.,

研究领域 (Research Area)

图片 (Image Data)



Western Blot analysis of rat kidney cells using CD300c Polyclonal Antibody. Antibody was diluted at 1:500. Secondary antibody was diluted at 1:20000



Western Blot analysis of RAT-kidney cells using CD300c Polyclonal Antibody diluted at 1: 500. Secondary antibody was diluted at 1:20000

注意事项 (Note)

For research use only .

产品名称: CD300c Rabbit Polyclonal Antibody
产品货号: APRab08341

