

产品名称: CD32-B/C Rabbit Polyclonal Antibody
产品货号: APRab08368



产品概述 (Summary)

产品名称 (Production Name)	CD32-B/C Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FCGR2B/FCGR2C FCGR2B; CD32; FCG2; IGFR2; Low affinity immunoglobulin gamma Fc region receptor II-b; IgG Fc receptor II-b; CDw32; Fc-gamma RII-b; Fc-gamma-RIIb;
别名 (Alternative Names)	FcRII-b; CD32; FCGR2C; CD32; FCG2; IGFR2;Low affinity immunoglobulin gamma Fc region receptor II-c; IgG Fc receptor II-c; CDw32; Fc-gamma RII-c; Fc-gamma-RIIc; FcRII-c; CD32
基因 ID (Gene ID)	2213.0
蛋白 ID (SwissProt ID)	P31994.The antiserum was produced against synthesized peptide derived from the Internal region of human FCGR2B/C. AA range:141-190

产品应用 (Application)

产品名称: CD32-B/C Rabbit Polyclonal Antibody
产品货号: APRab08368



稀释比 (Dilution Ratio) WB 1:500-1:2000, ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight) 35kDa


研究背景 (Background)

Low affinity immunoglobulin gamma Fc region receptor II-b/Low affinity immunoglobulin gamma Fc region receptor II-c encoded by FCGR2B is a low affinity receptor for the Fc region of immunoglobulin gamma complexes. The encoded protein is involved in the phagocytosis of immune complexes and in the regulation of antibody production by B-cells. Variations in this gene may increase susceptibility to systemic lupus erythematosus (SLE). Several transcript variants encoding different isoforms have been found for this gene. immune response,

研究领域 (Research Area)

B_Cell_Antigen;Fc gamma R-mediated phagocytosis;Systemic lupus erythematosus;

图片 (Image Data)



Western Blot analysis of K562 cells using CD32-B/C Polyclonal Antibody.. Secondary antibody was diluted at 1:20000

注意事项 (Note)

For research use only .