

产品名称: CD5 Rabbit Polyclonal Antibody
产品货号: APRab08411



产品概述 (Summary)

产品名称 (Production Name)	CD5 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD5
别名 (Alternative Names)	CD5; LEU1; T-cell surface glycoprotein CD5; Lymphocyte antigen T1/Leu-1; CD antigen CD5
基因 ID (Gene ID)	921.0
蛋白 ID (SwissProt ID)	P06127.The antiserum was produced against synthesized peptide derived from human CD5. AA range:421-470

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	60kDa

产品名称: CD5 Rabbit Polyclonal Antibody
产品货号: APRab08411




研究背景 (Background)

function: May act as a receptor in regulating T-cell proliferation. CD5 interacts with CD72/LYB-2., similarity: Contains 3 SRCR domains., function: May act as a receptor in regulating T-cell proliferation. CD5 interacts with CD72/LYB-2., similarity: Contains 3 SRCR domains.,


研究领域 (Research Area)

Hematopoietic cell lineage;

图片 (Image Data)



Western blot analysis of lysates from 293 cells, treated with PMA 125ng/ml 30', using CD5 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using CD5 Polyclonal Antibody

注意事项 (Note)

For research use only .