

产品名称: CD72 Rabbit Polyclonal Antibody
产品货号: APRab08447



产品概述 (Summary)

产品名称 (Production Name)	CD72 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CD72
别名 (Alternative Names)	CD72; B-cell differentiation antigen CD72; Lyb-2; CD72
基因 ID (Gene ID)	971.0
蛋白 ID (SwissProt ID)	P21854.The antiserum was produced against synthesized peptide derived from the Internal region of human CD72. AA range:170-220

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	40kDa

研究背景 (Background)

产品名称: CD72 Rabbit Polyclonal Antibody
产品货号: AP Rab08447




function: Plays a role in B-cell proliferation and differentiation. Associates with CD5., online information: CD72, similarity: Contains 1 C-type lectin domain., subunit: Homodimer; disulfide-linked., tissue specificity: Pre-B-cells and B-cells but not terminally differentiated plasma cells., function: Plays a role in B-cell proliferation and differentiation. Associates with CD5., online information: CD72, similarity: Contains 1 C-type lectin domain., subunit: Homodimer; disulfide-linked., tissue specificity: Pre-B-cells and B-cells but not terminally differentiated plasma cells.,


研究领域 (Research Area)

B_Cell_Antigen;

图片 (Image Data)




Western blot analysis of RAT-SPLEEN MOUSE-BRAIN MOUSE-SPLEEN lysis using CD72 antibody. Antibody was diluted at 1:2000. Secondary antibody was diluted at 1:20000




Western Blot analysis of various cells using Antibody diluted at 1:1000. Secondary antibody was diluted at 1:20000


产品名称: CD72 Rabbit Polyclonal Antibody
产品货号: APRab08447



Immunohistochemical analysis of paraffin-embedded Human Amygdala. 1, Antibody was diluted at 1:100 (4°,overnight) . 2, High-pressure and temperature EDTA, pH8.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200 (room temperature, 30min) .




Immunohistochemical analysis of paraffin-embedded Human Amygdala. 1, Antibody was diluted at 1:100 (4°,overnight) . 2, High-pressure and temperature EDTA, pH8.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200 (room temperature, 30min) .




Immunohistochemical analysis of paraffin-embedded Human Amygdala. 1, Antibody was diluted at 1:100 (4°,overnight) . 2, High-pressure and temperature EDTA, pH8.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200 (room temperature, 30min) .


产品名称: CD72 Rabbit Polyclonal Antibody
产品货号: APRab08447



Immunohistochemical analysis of paraffin-embedded Human spleen. 1, Antibody was diluted at 1:100 (4°,overnight) . 2, High-pressure and temperature EDTA, pH8.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200 (room temperature, 30min) .



Immunohistochemical analysis of paraffin-embedded Human spleen. 1, Antibody was diluted at 1:100 (4°,overnight) . 2, High-pressure and temperature EDTA, pH8.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200 (room temperature, 30min) .



Immunohistochemical analysis of paraffin-embedded Human spleen. 1, Antibody was diluted at 1:100 (4°,overnight) . 2, High-pressure and temperature EDTA, pH8.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200 (room temperature, 30min) .

产品名称: CD72 Rabbit Polyclonal Antibody
产品货号: APRab08447



注意事项 (Note)

For research use only .