

产品名称: CD85g Rabbit Polyclonal Antibody
产品货号: APRab08475



产品概述 (Summary)

产品名称 (Production Name)	CD85g Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	LILRA4 LILRA4; ILT7; Leukocyte immunoglobulin-like receptor subfamily A member 4;
别名 (Alternative Names)	CD85 antigen-like family member G; Immunoglobulin-like transcript 7; ILT-7; CD85g
基因 ID (Gene ID)	23547.0
蛋白 ID (SwissProt ID)	P59901.The antiserum was produced against synthesized peptide derived from the Internal region of human LILRA4. AA range:391-440

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	55kDa

产品名称: CD85g Rabbit Polyclonal Antibody
产品货号: APRab08475




研究背景 (Background)


This gene encodes an immunoglobulin-like cell surface protein that is expressed predominantly on plasmacytoid dendritic cells (PDCs) and modulates the function of these cells in the immune response. Expression of this gene is downregulated by interleukin 3 (IL3). This gene is one of a cluster of highly related genes located at chromosomal region 19q13.4. [provided by RefSeq, Jan 2015],function:May act as receptor for class I MHC antigens. Ligand binding leads to the activation of eosinophils and to the release of RNASE2, IL4 and leukotriene C4.,similarity:Contains 4 Ig-like C2-type (immunoglobulin-like) domains.,tissue specificity:Detected on eosinophils, neutrophils and monocytes.,

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of lysate from MCF7 cells, using LILRA4 Antibody.




Western Blot analysis of MCF7 cells using CD85g Polyclonal Antibody.. Secondary antibody was diluted at 1:20000

产品名称: CD85g Rabbit Polyclonal Antibody
产品货号: APRab08475



Immunohistochemical analysis of paraffin-embedded human-liver, antibody was diluted at 1:100



Immunohistochemical analysis of paraffin-embedded human-liver, antibody was diluted at 1:100

注意事项 (Note)

For research use only .