

产品名称: Cdc50B Rabbit Polyclonal Antibody
产品货号: APRab08525



产品概述 (Summary)

产品名称 (Production Name)	Cdc50B Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TMEM30B
别名 (Alternative Names)	TMEM30B; CDC50B; Cell cycle control protein 50B; Transmembrane protein 30B
基因 ID (Gene ID)	161291.0
蛋白 ID (SwissProt ID)	Q3MIR4.The antiserum was produced against synthesized peptide derived from human TMEM30B. AA range:91-140

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	38kDa

产品名称: Cdc50B Rabbit Polyclonal Antibody
产品货号: APRab08525




研究背景 (Background)

similarity:Belongs to the CDC50/LEM3 family.,similarity:Belongs to the CDC50/LEM3 family.,


研究领域 (Research Area)

Cell Biology


图片 (Image Data)



Western blot analysis of lysates from K562 cells, using TMEM30B Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HeLa cells using TMEM30B antibody.



Western Blot analysis of K562 cells using Cdc50B Polyclonal Antibody

产品名称: Cdc50B Rabbit Polyclonal Antibody
产品货号: APRab08525



注意事项 (Note)

For research use only .