

产品名称: Cdk11A/B Rabbit Polyclonal Antibody
产品货号: AP Rab08549



产品概述 (Summary)

产品名称 (Production Name)	Cdk11A/B Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CDK11A/CDK11B
别名 (Alternative Names)	CDK11A; CDC2L2; CDC2L3; PITSLREB; Cyclin-dependent kinase 11A; Cell division cycle 2-like protein kinase 2; Cell division protein kinase 11A; Galactosyltransferase-associated protein kinase p58/GTA; PITSLRE serine/threonine-protein kinase C
基因 ID (Gene ID)	728642/984
蛋白 ID (SwissProt ID)	Q9UQ88/P21127.Synthesized peptide derived from Cdk11A/B . at AA range: 190-270

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
----------------------	--------------------------------------

产品名称: Cdk11A/B Rabbit Polyclonal Antibody
产品货号: APRab08549




蛋白分子量 (Molecular Weight) 90kDa

研究背景 (Background)

This gene encodes a member of the serine/threonine protein kinase family. Members of this kinase family are known to be essential for eukaryotic cell cycle control. Due to a segmental duplication, this gene shares very high sequence identity with a neighboring gene. These two genes are frequently deleted or altered in neuroblastoma. The protein kinase encoded by this gene can be cleaved by caspases and may play a role in cell apoptosis. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2015], caution: The sequence shown here is derived from an Ensembl automatic analysis pipeline and should be considered as preliminary data.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of CDK11 Antibody. The lane on the right is blocked with the CDK11 peptide.

注意事项 (Note)

For research use only .